


ONLY ONE




Omni-Link

Family of Expansion Tubing
OMNI-LINK® & FLEX-LINK®



Resuscitating a patient can create compromising situations during transport. With that in mind, Mercury has developed the Omni-Link and Flex-Link family of expansion tubing to provide greater flexibility and longer reach capabilities. Each coupled with Mercury Medical's unique patient valve, the Omni-Link or Flex-Link does not add dead space to the patient circuit.



www.mercurymed.com

MercuryMedical®

Your Need ... Our Innovation®

Omni-Link®

Features

- Flexible corrugated tubing
- Available in 6" and 12" segments that can be interconnected to give you your desired length

Benefits

- Absorbs movement
- Provides greater versatility and longer reach capabilities
- Reduces arm fatigue and airway trauma
- Helps eliminate accidental disconnects and extubations

Flex-Link®

Features

- Flexible, expandable and collapsible tubing
- Available in collapsed 14" length, but can be lengthened up to 3' when expanded

Benefits

- Length desired for patient reach is defined by user
- Provides greater versatility and longer reach capabilities
- Reduces arm fatigue and potential for airway trauma
- Helps eliminate accidental disconnects and extubations
- Enables free movement of patient's head when intubated

CPR BAG



Synthetic Rubber

CPR-2 BAG



PVC

ORDERING INFORMATION

Omni-Link

Part #	Description	Packaged
#10-55360	Omni-Link (12" length)	50/Box
#10-55361	Omni-Link II (6" length)	25/Box

Flex-Link

Part #	Description	Packaged
#10-56024	Adult PVC CPR-2 Bag, Aero Reservoir, Mask, Manometer, Flex-Link and Bag Handle	6/Box
#10-56025	Adult PVC CPR-2 Bag, Aero Reservoir, Mask, Flex-Link and Bag Handle	6/Box
#10-56026	Adult PVC CPR-2 Bag, O ₂ Reservoir Bag, Mask, Manometer, Flex-Link and Bag Handle	6/Box
#10-56027	Adult PVC CPR-2 Bag, O ₂ Reservoir Bag, Mask, Flex-Link and Bag Handle	6/Box
#10-56034	Adult PVC CPR-2 Bag, O ₂ Reservoir Bag, Mask, Flex-Link, StatCO ₂	6/Box
#10-56040	Adult PVC CPR-2 Bag, No Bag Handle, O ₂ Reservoir, Flex-Link	6/Box
#10-56041	Adult PVC CPR-2 Bag, No Bag Handle, Mask, Aerosol Reservoir, Flex-Link	6/Box
#10-55079	Adult CPR Standard O ₂ Bag Reservoir, Intersurgical Quadralite™ Size 3 Mask, Flex-Link	6/Box
#10-56309	Adult PVC CPR-2 Bag, Expandable Reservoir, Mask, Manometer, Flex-Link and Bag Handle	12/Box
#10-56310	Adult PVC CPR-2 Bag, Expandable Reservoir, Mask, Flex-Link and Bag Handle	12/Box
#10-55362	Flex-Link (36" extended)	25/Box
#10-55363	Flex-Link (36" extended)	50/Box

11300 - 49th Street North • Clearwater, Florida 33762-4807 USA

Tel. (800) 237-6418 / (727) 573-0088 • Fax (800) 990-6375

www.mercurymed.com