

CAT® HOLDER

REV050918

Fits Securely and Easily on the Belt...

The C-A-T® Holder was specifically designed to allow personnel to place their C-A-T® Tourniquet on their vest or gear for rapid and easy access.

It is made of 1000D IR Signature Reduced Nylon with MOLLE/PALS-style connectors.

The C-A-T® Holder protects your tourniquet from the elements and has an easy-open elastic pull tab. *Accommodates both Gen 6 and Gen 7 C-A-T® Tourniquets (not included)*

Features:

- Tailored specifically for your C-A-T®, to protect it from exposure to the elements
- PALS/MOLLE-style connectors for versatile attachment to the duty vest, gear or belt (mounted horizontally or vertically) for rapid & easy access
- 1000D IR Signature Reduced Nylon
- Easy-open pull elastic tab with tourniquet ID patch – “TQ”
- NATO stock numbers available for COY and ODG. Contact Customer Service for details.

Dimensions:

- L 6.25 in. x W 1.75 in. x D 1.5 in.
- Weight: 1.26 oz

C-A-T® Holder

| Color | Item # | NSN# |
|-------|---------|------|
| BLK | 30-0057 | |
| COY | 30-0044 | |
| RED | 30-0128 | |
| ODG | 30-0053 | |
| MTC | 30-0059 | |



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®