

▲ COBALT® X



BLUE
NITRILE
NON-STERILE
AMBIDEXTROUS
FULLY TEXTURED
BEADED CUFF
FOR SINGLE USE ONLY
NON-LATEX

This thicker version of our popular Cobalt gloves offers the same exceptional fit and feel with the extra protection of a thicker glove. Fully textured to provide an exceptional grip in any condition.

• Extra thickness for additional protection





40% THICKER THAN OUR TRADITIONAL COBALT® GLOVE

PRODUCT CODE				
S	N211			
М	N212			
L	N213			
XL	N214			

PACKAGING			
100	Gloves/Box		
10	Boxes/Case		

GLOVE PROFILE		ASTM D 6319	COBALT®X
Length	mm (inches)	min. 220 (8.7)	240 (9.5)
Cuff Thickness	mm (mils)	n/a	0.10 (3.9)
Palm Thickness	mm (mils)	min. 0.05 (2.0)	0.13 (5.1)
Finger Thickness	mm (mils)	min. 0.05 (2.0)	0.16 (6.3)
Tensile Strength	Before Aging After Aging	min. 14 min. 14	27 14
Elasticity/Elongation	Before Aging After Aging	min. 500 min. 400	500 400
Pinhole Defect Rate*	(AQL)	max. 2.5	1.5

^{*} Pinhole Defect Rate is FDA minimum requirement.

















RECOMMENDED INDUSTRIES

