


SOFT STRETCHERS

BIG DADDY™



DMS **Soft Stretchers** are made with orange vinyl that is reinforced with black nylon strapping and vinyl handles.


- Rugged, yet comfortable vinyl handles added for easy gripping.
- Three sizes to choose from.

BIG DADDY™


60" X 80" – **BIG DADDY™** Soft Stretcher - weight load - 1,500 pounds
40" X 71" – **MEDIUM DADDY™** Soft Stretcher - weight load - 1,250 pounds
30" X 70" – **LITTLE DADDY™** Soft Stretcher - weight load - 1,000 pounds



CARRYING CASE INCLUDED WITH ALL SIZES



MEDIUM - 40" X 71"



SMALL - 30" X 70"

