

The Thomas Heating Pack is the largest and most complete IV fluid and medication warmer in the EMS industry.

Features of the pack include three interior color coded pockets with elastics designed to hold medications or other IV supplies, and two removable clear vinyl exterior pockets. The outside of the pack is constructed from 1000D water resistant fabric with protective foam and aluminum insulated polyethylene bubbles for maximum heat radiance and temperature maintenance. The exterior of the pack also includes reflective striping for improved visibility.

Additional Features:

- Adjustable back pack straps with chest connect for long distance rescue and treatment (Back pack straps can also be tucked away inside the pack if not needed)
- Holds up to 3 liters of IV fluids as well as all tubing accessories, medications, emergency blankets, blood pressure cuff and stethoscope, rope, and other rescue / treatment accessories
- Includes 1 fully insulated IV bag (pressure infused capable) with removable insulated tubing
- Fuse protected and easily removable ALUMINUM heating element
- Includes a redi-splint that serves as the frame of the pack and provides additional protection and insulation
- Pre-set 98 degree temperature maintained (+/- 1 degree) while plugged in
- Fuse protected and easily removable heating element
- 12 volt operation
- 1 year electronics warranty



Dimensions: 16" x 14" x 7"
 Heating pack comes complete with insulated IV bag and insulated IV tubing. Equipment not included



Insulated IV Pressure Infuser Bag



Insulated Tubing

Item	Description	Price
TT650	Thomas Heating Pack, Red	\$429.99
TT654	Thomas Heating Pack, Black	\$429.99

Item	Description	Price
TT521	Insulated IV Tubing - 72"	\$25.00
TT520	Insulated IV Pressure Infuser Bag	\$55.00
TT651	110 V Converter	\$55.00



Warranty: 1 year parts and labor.

Thomas EMS
 PO Box 651305
 Salt Lake City, Utah 84165
 Toll Free (800) 445-3640
 (801) 265-2516 • Fax (801) 268-9272
www.thomasems.com