




MedSource
Respiratory & Airway

ET Stylet



The MedSource ET Stylet is constructed of a malleable metal core that is coated in plastic, making it easy to shape. The stylet features a flexible, soft distal tip that reduces risk of patient trauma and a low-friction coating that makes it easy to withdraw from the endotracheal tube.

Packaged 20 per box.

Rev. 1

www.gomedsource.com

| Item # | Description | UOM |
|-----------|-----------------|-------------------------------|
| MS-ETS100 | ET Stylet, 6FR | 20 EA/BX, 10 BX/CS, 200 EA/CS |
| MS-ETS200 | ET Stylet, 10FR | 20 EA/BX, 10 BX/CS, 200 EA/CS |
| MS-ETS300 | ET Stylet, 12FR | 20 EA/BX, 10 BX/CS, 200 EA/CS |
| MS-ETS400 | ET Stylet, 14FR | 20 EA/BX, 10 BX/CS, 200 EA/CS |