



The Thomas Pediatric Pack can be used in conjunction with all existing pediatric systems.



Thomas EMS Fill Kit List for the Thomas Pedi Pack

3-5kg

| | |
|------------------------|-------------------------|
| 2.5 ET Tube | 24G x .75" IV Catheter |
| 3.0 ET Tube | 25G Butterfly Needle |
| 5 FR Suction Catheter | Miller #0 Disp. Blade |
| 6 FR Suction Catheter | IV Start Kit |
| 40mm Oral Airway | Alcohol Prep Pad |
| 50mm Oral Airway | Pediatric ECG Electrode |
| 10 FR Nasogastric Tube | |

6-9kg

| | |
|-------------------------|-------------------------|
| 2.5 ET Tube | 24G x .75" IV Catheter |
| 3.0 ET Tube | 23G Butterfly Needle |
| 3.5 ET Tube | 25G Butterfly Needle |
| 5 FR Suction Catheter | Miller #0 Disp. Blade |
| 6 FR Suction Catheter | Miller #1 Disp. Blade |
| 40mm Oral Airway | IV Start Kit |
| 50mm Oral Airway | Alcohol Prep Pad |
| 10 FR Nasogastric Tube | Pediatric ECG Electrode |
| 22G x 1.00" IV Catheter | |

10-11kg

| | |
|-------------------------|-------------------------|
| 4.0 ET Tube | 24G x .75" IV Catheter |
| 8 FR Suction Catheter | 21G Butterfly Needle |
| 10 FR Suction Catheter | 23G Butterfly Needle |
| 60mm Oral Airway | Miller #1 Disp. Blade |
| 10 FR Nasogastric Tube | Mac #1 Disp. Blade |
| 20G x 1.25" IV Catheter | IV Start Kit |
| 22G x 1.00" IV Catheter | Alcohol Prep Pad |
| | Pediatric ECG Electrode |

12-14kg

| | |
|-------------------------|-------------------------|
| 4.5 ET Tube | 21G Butterfly Needle |
| 10 FR Suction Catheter | 23G Butterfly Needle |
| 60mm Oral Airway | Miller #2 Disp. Blade |
| 10 FR Nasogastric Tube | Mac #2 Disp. Blade |
| 18G x 1.25" IV Catheter | IV Start Kit |
| 20G x 1.25" IV Catheter | Alcohol Prep Pad |
| 22G x 1.00" IV Catheter | Pediatric ECG Electrode |

15-18kg

| | |
|-------------------------|-------------------------|
| 5.0 ET Tube | 22G x 1.00" IV Catheter |
| 10 FR Suction Catheter | 21G Butterfly Needle |
| 60mm Oral Airway | 23G Butterfly Needle |
| 70mm Oral Airway | Miller #2 Disp. Blade |
| 10 FR Nasogastric Tube | Mac #2 Disp. Blade |
| 12 FR Nasogastric Tube | IV Start Kit |
| 18G x 1.25" IV Catheter | Alcohol Prep Pad |
| 20G x 1.25" IV Catheter | Pediatric ECG Electrode |

19-22kg

| | |
|-------------------------|-------------------------|
| 5.5 ET Tube | 21G Butterfly Needle |
| 10 FR Suction Catheter | 23G Butterfly Needle |
| 80mm Oral Airway | Miller #2 Disp. Blade |
| 12 FR Nasogastric Tube | Mac #3 Disp. Blade |
| 14 FR Nasogastric Tube | IV Start Kit |
| 18G x 1.25" IV Catheter | Alcohol Prep Pad |
| 20G x 1.25" IV Catheter | Pediatric ECG Electrode |

24-28kg

| | |
|-------------------------|-------------------------|
| 6.0 ET Tube | 21G Butterfly Needle |
| 10 FR Suction Catheter | Miller #2 Disp. Blade |
| 80 mm Oral Airway | Miller #3 Disp. Blade |
| 14 FR Nasogastric Tube | Mac #3 Disp. Blade |
| 16 FR Nasogastric Tube | Mac #4 Disp. Blade |
| 18G x 1.25" IV Catheter | IV Start Kit |
| 20G x 1.25" IV Catheter | Alcohol Prep Pad |
| 19G Butterfly Needle | Pediatric ECG Electrode |

30-36 kg

| | |
|-------------------------|-------------------------|
| 6.5 ET Tube | 19G Butterfly Needle |
| 12 FR Suction Catheter | 21G Butterfly Needle |
| 90 mm Oral Airway | Miller #3 Disp. Blade |
| 18 FR Nasogastric Tube | Mac #4 Disp. Blade |
| 16G x 1.25" IV Catheter | IV Start Kit |
| 18G x 1.25" IV Catheter | Alcohol Prep Pad |
| 20G x 1.25" IV Catheter | Pediatric ECG Electrode |

Non-weight specific equipment

500 ml Sodium Chloride
3-way Stopcock
10 Drop – IV Set
60 Drop – IV Set
1" Tape
Infant – BP Cuff
Child – BP Cuff
Infant – Cervical Collar

Pediatric – Cervical Collar
Pediatric - Stethoscope
Infant – Nasal Cannula
Pediatric – Nasal Cannula
Pediatric – CO2 Detector
Pediatric – Magill Forceps
2 oz. Bulb Syringe
15G – I.O. Needle

18G – I.O. Needle
Sharps Shuttle
6" IV Arm Board
9" IV Arm Board
Small - Laryngoscope Handle
Medium - Laryngoscope Handle
Infant - BVM
Child - BVM

Neonate - Disposable Mask
Infant - Disposable Mask
Toddler - Disposable Mask
Small - Disposable Mask
Infant - 02 Non-rebreather Mask
Pediatric – 02 Non-rebreather Mask