


# JELCO®

## I.V. CATHETER

THE WORLDWIDE "GOLD STANDARD" IN I.V. THERAPY



### Sharp, beveled needle tip delivers consistent, easy insertions

- Thinwall catheter design provides higher flow rates allowing for the use of a smaller gauge device
- FEP polymer catheter and J-point needle design allow clinician to "feel" vein entry
- Excellent trim accuracy allows smoother vein access
- JELCO®-W I.V. Catheter with "WINGS," facilitates hub securement



### Trusted by healthcare professionals for over 30 years

- Pediatric/neonatal communities appreciate JELCO® I.V. Catheters for their high flow rates and improved securement for the most active patients
- Anesthesiologists rely on JELCO® I.V. Catheters for their consistent flash, high flow rates, and ease of use for arterial lines and gaining I.V. access
- A high-performance product plus the confidence in knowing it comes from Vascular Access


smiths

Smiths Medical

# JELCO®

## I.V. CATHETER

THE WORLDWIDE "GOLD STANDARD" IN I.V. THERAPY



Catheter gauge and Length	JELCO® I.V. Catheter <i>Radiopaque FEP Polymer, Straight Hub</i>	JELCO® I.V. Catheter <i>Striped, Radiopaque FEP Polymer, Straight Hub</i>	JELCO®-W I.V. Catheter <i>Radiopaque FEP Polymer, Winged Hub</i>
14G x 1 1/4"	4048	—	4078
14G x 2"	4088	4098	—
16G x 1 1/4"	4042	—	4072
16G x 2"	4082	4092	—
18G x 1 1/4"	4055	4065	4075
18G x 1 3/4"	4054	4064	—
20G x 1"	4057	4067	4077
20G x 1 1/4"	4056	4066	4076
20G x 1 3/4"	4059	4069	—
22G x 1"	4050	4060	4070
24G x 5/8"	—	—	4073
24G x 3/4"	4053	4063	—

50 per box, 4 boxes per case

*I.V. catheters may vary greatly with regard to materials, performance, and safety features. This guide is designed primarily for cross-referencing I.V. catheter gauges and lengths. Please contact your local Smiths Medical representative for assistance in cross-referencing to the I.V. catheter that best fits your clinical needs.*

**Smiths Medical**  
2231 Rutherford Road  
Carlsbad, CA 92008  
USA

T: 800 848 1757  
F: 800 621 2590

©2005 Smiths Medical family of companies. All rights reserved. Medex and Jelco are trademarks of the Smiths Medical family of companies. The ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.

smiths

Smiths Medical

VA405 2005.5