



WE CARE LIKE FAMILY™



Tracheostomy Tube Holders

Keeps tracheostomy tube securely in place

- Soft strap features an adjustable hook-and-loop closure
- One size fits all
- Soft compliant foam laminate with elastic adjustment increases patient comfort and accommodates coughing and edema while still providing a snug fit
- Soft twill ribbon between collar and tube helps prevent sliding and patient irritation, and also helps prevent neckband from becoming soiled
- Designed as an aid in preventing accidental ventilator disconnections at the trach site
- Also stabilizes the trach tube
- Vent tie; anti-disconnect
- Not made with natural rubber latex



4620



4621



Reorder No.	Description	Quantity
4620	Trach Tube Holder - Adult (19" length)	1/polybag 5/10/Cs
4621	Trach Tube Holder - Pediatric (9" length)	1/polybag 5/10/Cs

HCPCS Code-A7526



See website for our Tracheostomy Care Kits (Reorder No. 4600-4602)

dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

© 2016 Dynarex Corporation



Find us on:

