

Product Name	qubeSERIES Monitor Simulation
Training Curricula	American Heart Association (AHA) European Resuscitation Council (ERC) National Association of Emergency Technicians (NAEMT) International Trauma Life Support (ITLS) Advanced Trauma Life Support (ATLS) National Association of State EMS Officials (NASEMSO) National Registry of Emergency Medical Technicians (NREMT)
Product Name/ Order Codes	<p>Patient Monitors/Defibrillators Simulators: qube 15 (LIFEPAK 15-STRYKER) SKU#: SQqube15 qube X (X-Series - ZOLL) SKU#: SQqubeX qube3 (Corpuls3 simulation) SKU#: SQqube3 qube3T (Corpuls3T simulation) SKU#: SQqube3T qube7 (DEFIGARD Touch7 - Schiller) SKU#: SQqube7 qubeZERO (Generic) SKU#: SQqubeZero qube20e (LIFEPAK 20e in hospital simulator) SKU#: SQqube20e qubeR (Zoll R-Series in hospital simulator) SKU#: SQqubeR qubeMC (MEDUCORE - WEINMANN) SKU#: SQqubeMC</p> <p>Ventilator: qubeMS2 (MEDUMAT Standard2 - WEINMANN) SKU#: SQqubeMS2 qubeVENT (Generic) SKU#: SQqubeVENT</p> <p>AV Capture & Debriefing: qubeAVgo SKU#: SQqubeAVgo qubeAVpro SKU#: SQqubeAVpro</p>
Note	Choice of hardware sets packaged with the product for order fulfilment will vary based on restrictions of Apple products used in the configuration.

Disclaimer: These devices are intended to be used as an instructional aide in the teaching and demonstrating of training. Actual products may vary slightly from photos. Nasco Healthcare reserves the right to change product colors, materials, supplies, or functions as needed.



Nasco Healthcare Inc.
 16 Simulaid Drive, Saugerties, NY, 12477 USA
 Phone: 1-833-NASCOHC (627-2642) Email: info@nascohealthcare.com
 www.nascohealthcare.com



Be
READY

Monitor Simulation: SKILLQUBE



MONITOR TRAINING, AS REALISTIC AS POSSIBLE.



qubeSERIES MONITOR INTERFACE

Powerful simulation systems for your patient monitor/defibrillators. SKILLQUBE offers a wide array of Patient Monitor/Defibrillator options designed to support different named medical device manufacturers available in the market (Stryker, Zoll, Weinmann and Schiller); a generic version is also available for purchase.

- Full functional ECG-monitor interface
- CPR-mode, defibrillation, pacemaker
- Checklists and standard procedures for documentation
- Seamless integration with qubeAV mobile debriefing system
- Browser-based scenario editor, image and video library
- Control of vital parameters by instructor
- Auscultation functionality, blood glucometer, thermometer with an assessment app



qubeZERO



qube³



qube^{3T}



qube⁷



qube¹⁵



qube^X



qube^{AED}



qube^{20e}



qube^R



qube^{MC}

VENTILATION SYSTEM

The MEDUMAT Standard2 is used in emergency services, air rescue or in the military's medical service.



- Mode for resuscitation (CPR)
- Simulation of flow measurement and CO₂ measurement
- Monitoring of IPPV, CPAP and oxygen administration
- Instructor controls for iPad
- Preset and custom scenarios



qube^{MS2}



qube^{VENT}



qubeAV DEBRIEFING SYSTEM

Debriefing is one of the most important parts of simulation training. With our qubeAV debriefing systems, you will be able to support students during training and re-run the debriefing after the training.



Mobile debriefing-system with audio-video-recording for simulation training



Due to the mobility of the system, you are totally flexible, and you can choose if you like to train in the classroom, on the ward, or in the ambulance



Change of perspective: flexible change between up to four cameras during training



Seamless integration with all patient simulation systems



Recording the training allows analysis of the training scenario from all angles



Multiple marker-sets defined by you, controlled via the marker keyboard, track pad or iPad